

[Home Inspections, Compliance, Enforcement, and Criminal Investigations](#) [Inspections](#) [Inspection Guides](#)

Inspections, Compliance, Enforcement, and Criminal Investigations

Measurement of Relative Humidity in the ETO Process

[\[Previous Chapter¹\]](#) [\[Table of Contents²\]](#) [\[Next Chapter³\]](#)

**DEPT. OF HEALTH, EDUCATION,
AND
WELFARE PUBLIC HEALTH
SERVICE
FOOD AND DRUG ADMINISTRATION
*ORA/ORO/DEIO/IB***

**Date: 4/30/87 Number: 47
Related Program Areas:
Drugs and Medical Devices**

ITG SUBJECT: MEASUREMENT OF RELATIVE HUMIDITY IN THE EO PROCESS

The importance and criticality of the moisture content or relative humidity (RH) of the environment within the enclosed vessel of the Ethylene Oxide (EO) sterilization process has been well established by the drug and device industries. Appropriate instrumentation systems are necessary to sense, control, and record this parameter. Of these, the key element of the system and the most difficult to establish and maintain in terms of accuracy is the moisture sensor. This ITG will briefly describe basic types of sensors together with some limitations.

TYPES OF SENSORS/INSTRUMENTS

The simplest form of RH measurement and one with which most are familiar, is the traditional Mechanical Hygrometer. It is based on the principle that certain materials will expand or contract as a function of the amount of moisture absorbed. Materials such as paper, wood, bone, plant leaves, textiles, animal tissue and plastics have been historically used. The most common materials, however, have been human or horse hair. Direct readout is usually obtained on a dial indicator interconnected to the material via mechanical linkage. It is also obvious that such an instrument would be difficult for remote reading in the enclosed space of a sterilization chamber and for all practicality could not be used. It is also obvious that precision would not be a relied-upon quality. Response is slow for both increasing and decreasing humidity and temperature compensation cannot be made (calibration applies only to a specific temperature). Also, the aging process of the material is not uniform, making calibration extremely difficult. It is doubtful that the investigator will ever encounter the use of such an instrument to measure the RH in an EO chamber.

Electric Hygrometer - This type utilizes a sensor that is a hygroscopic film. An instantaneous change of electrical resistance or capacitance can be measured as the result of small changes in absorbed moisture by the film. The hygroscopic film is sometimes replaced by a wire grid wrapped around a substrate. The wire grid is usually coated with a hygroscopic salt solution (lithium chloride is common) which absorbs water vapor. Resistance of the wire grid corresponds to a calibrated RH value. A calibration curve is usually furnished for each unit by the manufacturer. In physical appearance for EO chamber use, the sensor most likely will resemble a metallic or plastic cylinder approximately 3/4 inch diameter by 2 inch long with perforations on the side and electrical contact pins or connector on an end. The response to changing RH is rapid and accuracy tolerances can be as low as $\pm 1.5\%$ RH. This type of sensor usually is not reliable after multiple exposures to EO. A maintenance program to clean, rejuvenate, and recalibrate periodically is necessary for the user. The investigator should also be alerted to the fact that humidity sensitive materials are also temperature dependent. This means that the sensor must be calibrated over the entire temperature range of operation.

Dew-Point Hygrometer - This is commonly known as the chilled mirror technique wherein the optical qualities of a mirror are measured as it is cooled by the sampled gas to cause condensation. The amount of condensation is related to the dew-point temperature which together with a dry bulb temperature and a standard chart can be interpolated to reveal RH. This method can be quite accurate ($\pm 2\%$ RH). One manufacturer utilizes microprocessor technology to offer simultaneous measurement of not only dew-point but temperature and RH as well, with a digital readout display. Accuracy in determining RH has been claimed as low as $\pm 0.5\%$. Maintenance required for this type of sensor would be primarily the cleaning of the mirror and sampled gas line.

Gas Liquid Chromatography (GLC) - This technique utilizes sampling of the gas mixtures in the contained vessel. The area under the curves generated by the GLC are measured manually or automatically with a microprocessor to determine concentrations of sterilant gases, air, and moisture contents (RH). This is probably the most accurate means to measure RH. The critical factor is proper maintenance of the instrument and periodic calibration. The sampling ports must be clear and free of any debris and the sampling lines must be heated and insulated to prevent condensation. Not everyone understands the principles of GC; therefore, a qualified person, usually a chemist should be responsible for adjustments or maintenance of the equipment. This method also lends itself quite well to process control systems.

CONCLUSION

Other forms of sensors and instruments are available for measuring RH. Most are variations of the basic types described but they offer much less in practicality for the EO sterilizer. Some firms choose to use an empirical method for determining RH. Based upon the laws of thermodynamics for an ideal gas, RH, which is a ratio, can be expressed as a percentage of the mole fraction of water vapor in a space (chamber) to the mole fraction of water vapor in the space at saturation at a particular temperature. In the case of the EO chamber, this could be expressed in pressures for practical engineering applications by the equation:

$$\% \text{ RH} = \frac{\Delta P}{P_{sv}} \times 100$$

P sv (Saturation vapor pressure) psia

Where, Δ P is the change in absolute pressure of the chamber resulting from steam injection after the pre-evacuation portion of the cycle and P sv is the saturation vapor pressure taken from published tables of saturated steam and temperature. (See the attached excerpt of such a table). As an example, steam is injected into the EO chamber to increase the pressure by 3.0 in. Hg. (1.47 psia). Sterilization temperature is 130 F and from the attached table of saturation vapor pressure at a temperature, P sv = 2.2230 psia. The equation then becomes:

$$\text{RH} = \frac{1.47}{2.2230}$$

$$\text{-----} \times 100 = 66\%$$

$$2.2230$$

The use of the empirical method is adequate providing that all instrumentation to measure both pressures and temperature are adequately calibrated and maintained. Validation/qualification of the process should verify the accuracy of the calculations. Some firms utilize this method as a means for back-up or for cross-checking to determine faulty sensors. The investigator is invited to also utilize this method for a rapid check of the RH parameter of the EO process during the Establishment Inspection.

References

1. The American Society of Mechanical Engineers; ASME Steam Tables; 1967.
2. Robertson, J.H., Townsend, M.W., Allen, P.M., Devisser, A., and Enzinger, R.M.; Validation of Ethylene Oxide Sterilization Cycles; PDA Spring Meeting; Chicago, IL; March 25, 1977.
3. Quin, F.C.; Humidity/Moisture Considerations; American Instrument Co., Silver Spring, MD; Reprint No. 490.

i Attachment to ITG #47

Saturated Steam: Temperature Table
 Chamber
 Temp. Psv

Abs Press. Specific Volume Enthalpy Entropy

Temp	Lb per	Sat.	Sat.	Sat.	Sat.	Sat.	Sat.	Sat.	Sat.	Sat.	Temp
Fahr	Sq In.	Liquid	Evap	Vapor	Liquid	Evap	Vapor	Liquid	Evap	Vapor	Fahr
t	p	v f	v fg	v g	h f	h fg	h g	S f	S fg	S g	t
32.0*	0.08859	0.016022	3304.7	3304.7	-0.0179	1075.5	1075.5	0.0000	2.1873	2.1873	32.0*
34.0	0.09600	0.016021	3061.9	3061.9	1.996	1074.4	1076.4	0.0041	2.1762	2.1802	34.0
36.0	0.10395	0.016020	2839.0	2839.0	4.008	1073.2	1077.2	0.0081	2.1651	2.1732	36.0
38.0	0.11249	0.016019	2634.1	2634.2	6.018	1072.1	1078.1	0.0122	2.1541	2.1663	38.0
40.0	0.12163	0.016019	2445.8	2445.8	8.027	1071.0	1079.0	0.0162	2.1432	2.1594	40.0
42.0	0.13143	0.016019	2272.4	2272.4	10.035	1069.8	1079.9	0.0202	2.1325	2.1527	42.0
44.0	0.14192	0.016019	2112.8	2112.8	12.041	1068.7	1080.7	0.0242	2.1217	2.1459	44.0
46.0	0.15314	0.016020	1965.7	1965.7	14.047	1067.6	1081.6	0.0282	2.1111	2.1393	46.0
48.0	0.16514	0.016021	1830.0	1830.0	16.051	1066.4	1082.5	0.0321	2.1006	2.1327	48.0
50.0	0.17796	0.016023	1704.8	1704.8	18.054	1065.3	1083.4	0.0361	2.0901	2.1262	50.0
52.0	0.19165	0.016024	1589.2	1589.2	20.057	1064.2	1084.2	0.0400	2.0798	2.1197	52.0
54.0	0.20625	0.016026	1482.4	1482.4	22.058	1063.1	1085.1	0.0439	2.0695	2.1134	54.0
56.0	0.22183	0.016028	1383.6	1383.6	24.059	1061.9	1086.0	0.0478	2.0593	2.1070	56.0
58.0	0.23843	0.016031	1292.2	1292.2	26.060	1060.8	1086.9	0.0516	2.0491	2.1008	58.0
60.0	0.25611	0.016033	1207.6	1207.6	28.060	1059.7	1087.7	0.0555	2.0391	2.0946	60.0
62.0	0.27494	0.016036	1129.2	1129.2	30.059	1058.5	1088.6	0.0593	2.0291	2.0885	62.0
64.0	0.29497	0.016039	1056.5	1056.5	32.058	1057.4	1089.5	0.0632	2.0192	2.0824	64.0
66.0	0.31626	0.016043	989.0	989.1	34.056	1056.3	1090.4	0.0670	2.0094	2.0764	66.0
68.0	0.33889	0.016046	926.5	926.5	36.054	1055.2	1091.2	0.0708	1.9996	2.0704	68.0
70.0	0.36292	0.016050	868.3	868.4	38.052	1054.0	1092.1	0.0745	1.9900	2.0645	70.0
72.0	0.38844	0.016054	814.3	814.3	40.049	1052.9	1093.0	0.0783	1.9804	2.0587	72.0
74.0	0.41550	0.016058	764.1	764.1	42.046	1051.8	1093.8	0.0821	1.9708	2.0529	74.0
76.0	0.44420	0.016063	717.4	717.4	44.043	1050.7	1094.7	0.0858	1.9614	2.0472	76.0
78.0	0.47461	0.016067	673.8	673.9	46.040	1049.5	1095.6	0.0895	1.9520	2.0415	78.0

80.0	0.50683	0.016072	633.3	633.3	48.037	1048.4	1096.4	0.0932	1.9426	2.0959	80.0
82.0	0.54093	0.016077	595.5	595.5	50.033	1047.3	1097.3	0.0969	1.9334	2.0303	82.0
84.0	0.57702	0.016082	560.3	560.3	52.029	1046.1	1098.2	0.1006	1.9242	2.0248	84.0
86.0	0.61518	0.016087	227.5	527.5	54.026	1045.0	1099.0	0.1043	1.9151	2.0193	86.0
88.0	0.65551	0.016093	496.8	496.8	56.022	1043.9	1099.9	0.1079	1.9060	2.0139	88.0
90.0	0.69813	0.016099	468.1	468.1	58.018	1042.7	1100.8	0.1115	1.8970	2.0086	90.0
92.0	0.74313	0.016105	441.3	441.3	60.014	1041.6	1101.6	0.1152	1.8881	2.0033	92.0
94.0	0.79062	0.016111	416.3	416.3	62.010	1040.5	1102.5	0.1188	1.8792	1.9980	94.0
96.0	0.84072	0.016117	392.8	392.9	64.006	1039.3	1103.3	0.1224	1.8704	1.9928	96.0
98.0	0.89356	0.016123	370.9	370.9	66.003	1038.2	1104.2	0.1260	1.8617	1.9876	98.0
100.0	0.94924	0.016130	350.4	350.4	67.999	1037.1	1105.1	0.1295	1.8530	1.9825	100.0
102.0	1.00789	0.016137	331.1	331.1	69.995	1035.9	1105.9	0.1331	1.8444	1.9775	102.0
104.0	1.06965	0.016144	313.1	313.1	71.992	1034.8	1106.8	0.1366	1.8358	1.9725	104.0
106.0	1.1347	0.016151	296.16	296.18	73.99	1033.6	1107.6	0.1402	1.8273	1.9675	106.0
108.0	1.2030	0.016158	280.28	280.30	75.98	1032.5	1108.5	0.1437	1.8188	1.9626	108.0
110.0	1.2750	0.016165	265.37	265.39	77.98	1031.4	1109.3	0.1472	1.8105	1.9577	110.0
112.0	1.3505	0.016173	251.37	251.38	79.98	1030.2	1110.2	0.1507	1.8021	1.9528	112.0
114.0	1.4299	0.016180	238.21	238.22	81.97	1029.1	1111.0	0.1542	1.7938	1.9480	114.0
116.0	1.5133	0.016188	225.84	225.85	83.97	1027.9	1111.9	0.1577	1.7856	1.9433	116.0
118.0	1.6009	0.016196	214.20	214.21	85.97	1026.8	1112.7	0.1611	1.7774	1.9386	118.0
120.0	1.6927	0.016204	203.25	203.26	87.97	1025.6	1113.6	0.1646	1.7693	1.9339	120.0
122.0	1.7891	0.016213	192.94	192.95	89.96	1024.5	1114.4	0.1680	1.7613	1.9293	122.0
124.0	1.8901	0.016221	183.23	183.24	91.96	1023.3	1115.3	0.1715	1.7533	1.9247	124.0
126.0	1.9959	0.016229	174.08	174.09	93.96	1022.2	1116.1	0.1749	1.7453	1.9202	126.0
128.0	2.1068	0.016238	165.45	165.47	95.96	1021.0	1117.0	0.1783	1.7374	1.9157	128.0
130.0	2.2230	0.016247	157.32	157.33	97.96	1019.8	1117.8	0.1817	1.7295	1.9112	130.0
132.0	2.3445	0.016256	149.64	149.66	99.95	1018.7	1118.6	0.1851	1.7217	1.9068	132.0
134.0	2.4717	0.016265	142.40	142.41	101.95	1017.5	1119.5	0.1884	1.7140	1.9024	134.0
136.0	2.6047	0.016274	135.55	135.57	103.95	1016.4	1120.3	0.1918	1.7063	1.8980	136.0

138.0 2.7438 0.016284 129.09 129.11 105.95 1015.2 1121.1 0.1951 1.6986 1.8937 138.0

140.0 2.8892 0.016293 122.98 123.00 107.95 1014.0 1122.0 0.1985 1.6910 1.8895 140.0

142.0 3.0411 0.016303 117.21 117.22 109.95 1012.9 1122.8 0.2018 1.6534 1.8852 142.0

144.0 3.1997 0.016312 111.74 111.76 111.95 1011.7 1123.6 0.2051 1.6759 1.8810 144.0

146.0 3.3653 0.016322 106.58 106.59 113.95 1010.5 1124.5 0.2084 1.6684 1.8769 146.0

148.0 3.5381 0.016332 101.68 101.70 115.95 1009.3 1125.3 0.2117 1.6610 1.8727 148.0

150.0 3.7184 0.016343 97.05 97.07 117.95 1008.2 1126.1 0.2150 1.6536 1.8686 150.0

152.0 3.9065 0.016353 92.66 92.68 119.95 1007.0 1126.9 0.2183 1.6463 1.8646 152.0

154.0 4.1025 0.016363 88.50 88.52 121.95 1005.8 1127.7 0.2216 1.6390 1.8606 154.0

156.0 4.3068 0.016374 84.56 84.57 123.95 1004.6 1128.6 0.2248 1.6318 1.8566 156.0

158.0 4.5197 0.016384 80.82 80.83 125.96 1003.4 1129.4 0.2281 1.6245 1.8526 158.0

160.0 4.7414 0.016395 77.27 77.29 127.96 1002.2 1130.2 0.2313 1.6174 1.8487 160.0

162.0 4.9722 0.016406 73.90 73.92 129.96 1001.0 1131.0 0.2345 1.6103 1.8448 162.0

164.0 5.2124 0.016417 70.70 70.72 131.96 999.8 1131.8 0.2377 1.6032 1.8409 164.0

166.0 5.4623 0.016428 67.67 67.68 133.97 998.6 1132.6 0.2409 1.5961 1.8371 166.0

168.0 5.7223 0.016440 64.78 64.80 135.97 997.4 1133.4 0.2441 1.5892 1.8333 168.0

170.0 5.9926 0.016451 62.04 62.06 137.97 996.2 1134.2 0.2473 1.5822 1.8295 170.0

172.0 6.2736 0.016463 59.43 59.45 139.98 995.0 1135.0 0.2505 1.5753 1.8258 172.0

174.0 6.5656 0.016474 56.95 56.97 141.98 993.8 1135.8 0.2537 1.5684 1.8221 174.0

176.0 6.8690 0.016486 54.59 54.61 143.99 992.6 1136.6 0.2568 1.5616 1.8184 176.0

178.0 7.1840 0.016498 52.35 52.36 145.99 991.4 1137.4 0.2600 1.5548 1.8147 178.0

i

[\[Previous Chapter⁴\]](#) [\[Table of Contents⁵\]](#) [\[Next Chapter⁶\]](#)

Page Last Updated: 02/17/2009

Note: If you need help accessing information in different file formats, see [Instructions for Downloading Viewers and Players](#).

[Accessibility Contact](#) [FDA Careers](#) [FDA Basics](#) [FOIA No Fear Act](#) [Site Map](#) [Transparency](#) [Website Policies](#)

U.S. Food and Drug Administration
10903 New Hampshire Avenue
Silver Spring, MD 20993

Ph. 1-888-INFO-FDA (1-888-463-6332)

[Email FDA](#)



[For Government](#) [For Press](#)

[Combination Products](#) [Advisory Committees](#) [Science & Research](#) [Regulatory Information](#) [Safety Emergency](#)
[Preparedness](#) [International Programs](#) [News & Events](#) [Training and Continuing Education](#) [Inspections/Compliance](#)
[State & Local Officials](#) [Consumers](#) [Industry](#) [Health Professionals](#) [FDA Archive](#)



U.S. Department of **Health & Human Services**

Links on this page:

1. <http://wcms.fda.gov/FDAgov/ICECI/Inspections/InspectionGuides/InspectionTechnicalGuides/UCM072925.htm>
2. <http://wcms.fda.gov/FDAgov/ICECI/Inspections/InspectionGuides/InspectionTechnicalGuides/default.htm>
3. [ssLINK/UCM081215.htm](#)
4. <http://wcms.fda.gov/FDAgov/ICECI/Inspections/InspectionGuides/InspectionTechnicalGuides/UCM072925.htm>
5. <http://wcms.fda.gov/FDAgov/ICECI/Inspections/InspectionGuides/InspectionTechnicalGuides/default.htm>
6. [ssLINK/UCM081215.htm](#)